

IMPORTANT! Please Read

Conventional serological pipetting

When you first receive your ali-Q™ 2, press the purple button to wake it out of shipping mode.



Note:

1. The motor will activate at seemingly random times when pressing the aspirate or dispense buttons. Don't worry. It's normal!
2. ali-Q 2 will go into sleep mode after 3 minutes of non-use but will wake with any movement of the unit.

For more information on ali-Q's intelligent measuring system, see Operating Manual at www.vistalab.com/support

Aliquoting/repeat pipetting


speed dial




ali-Q 2 VS model only:

Set desired aliquoting speed using speed dial

aliquot dials



All ali-Q 2 models:

1. Set aliquot volume on aliquot dial.
2. Aspirate adequate volume for total number of desired aliquots.
3. Press and HOLD the purple  aliquot button until the aliquot is dispensed. You will hear and feel the motor come on. Repeat press and hold for each aliquot. For maximum precision, discard first aliquot.

aliquot button



Note: if you don't hold the button long enough for the full aliquot volume, you will hear an "alert" sound and red light will flash.



Description

Description	VistaLab Catalog #
ali-Q 2 pipet controller	2100-0005
ali-Q 2 LS pipet controller	2100-0003
ali-Q 2 VS 5.0 mL pipet controller	2100-1005
ali-Q 2 VS 10.0 mL pipet controller	2100-1010

Additional box contents: spare	9020-4003
0.22µm membrane filter power	9020-9005
supply w/outlet adapters wall mount	9020-9008
calibration certificate	

**Register your ali-Q 2
online at
www.vistalab.com
for 1 year Warranty**

Go to www.vistalab.com/support
for full Operating Manual



**Complete your optimal liquid handling system with
Wobble-not™ serological pipets**

Standard Length Catalog

5 mL (200/cs)	4090-0005
10 mL (200/cs)	4090-0010
25 mL (200/cs)	4090-0025
50 mL (100/cs)	4090-0050
100 mL (50/cs)	4090-0100

Shorty Length

10 mL (200/cs)	5090-0010
25 mL (100/cs)	5090-0025



Think Differently. Feel Better.

2 Geneva Road Brewster, NY 10509 USA
(888) 652-6520; (914) 244-6226
customer.service@vistalab.com

©2020 VistaLab Technologies, Inc. All rights reserved.

L080-0052-001 RevC