

MLA Adjustable Volume Models	VistaLab Catalog #	Model	Color	Volume Range	Dispensing Increments	Accuracy*	Precision*
	1132	D-Tipper 50/100/200µL selectable pipette	Silver	50–200µL		±2.5% at 50µL ±1.5% at 100µL ±1.3% at 200µL	0.80% at 50µL 0.54% at 100µL 0.30% at 200µL
	1061	Macro 1–5mL adjustable pipette	Blue	1-5mL		±1.0% at 5mL ±1.0% at 3mL ±1.0% at 1mL	0.60% at 5mL 0.50% at 3mL 0.50% at 1mL
	1134	Macro 2,3,4,5 mL selectable pipette	Blue	2-5mL		±1.0%	0.50%
	1136	Digital pipette 10–50µL	Green	10–50µL	0.5µL	±5% at 10µL ±3% at 30µL ±2% at 50µL	2.6% at 10µL 0.93% at 30µL 0.80% at 50µL
	1137	Digital pipette 50–200µL	Red	50–200µL	1.0µL	±2% at 50µL ±1% at 125µL ±1% at 200µL	0.80% at 50µL 0.43% at 125µL 0.50% at 200µL
	1138	Digital pipette 200–1000µL	Gray	200–1000µL	5.0µL	±1% at 200µL ±1% at 600µL ±1% at 1000µL	0.50% at 200µL 1.00% at 600µL 0.60% at 1000µL

