

FREE!

Wobble-notTM
Low Insertion Force Serological Pipet

FREE!

Buy an Ovation **Macro10** electronic pipette
Get **2 FREE** cases of 10 mL Wobble-not serological pipettes



Macro10

10 mL macro volume
Electronic Pipette

VistaLab Catalog # 1065-10mL

Ovation® Macro10

High precision alternative to portable pipette controllers

- Uses Ovation 10 mL macro tips or serological pipets
- Fully electronic operation
- 2 Custom programmable volume pre-sets
- Programmable functions for:
 - multi-dispensing
 - serial dilutions
 - mixing
 - reverse pipetting
 - pipetting speed
- Power supply included



Demo Program

Ask about our
No Obligation
2 week trial

or see

www.vistalab.com/request-a-demo

To find a local distributor, go to www.vistalab.com/distributors

2 cases of Wobble-not 10 mL serological pipets

FREE

with Macro10 purchase

**Hassle-Free
Pipetting!**

**No Drip
Design!**



2-tier step design

- = 2 zones of contact in your pipet controller
- = Less force to insert
- = Less "wobble" = Max Stability
- = Better seal = Less dripping = Better Accuracy

Find your VistaLab Regional Manager



- Keith Brodnick, Canada
(847) 521-0325, kbrodnick@vistalab.com
- Gary Steinke, Northeast
(845) 803-1012, gsteinke@vistalab.com
- Drew Hamilton, Southeast
(845) 661-2342, dhamilton@vistalab.com
- Marie Ferguson, Central
(914) 244-4048, mferguson@vistalab.com
- Joe Tracy, West
(619) 895-6865, jtracy@vistalab.com

Cannot be combined with any other offer. Macro10 must be purchased by 12/31/16 and redemption request for free Wobble-nots must be received by 1/16/17. To receive free Wobble-not, contact your VistaLab Regional Manager or customer.service@vistalab.com. You must provide proof of Macro10 purchase. Reference promo code macro10wn.

Think Differently. Feel Better.