

**Q4 Special
Promo**



Trade-in your old pipettes
and give up the pain!
Then...

Buy an Ovation Multichannel and Get an Ovation M FREE or choose \$250 in FREE product*

Multichannel pipette must be purchased between 10/1/18 - 12/31/18.
Proof of purchase required. Free product must be claimed by 1/15/19.
Currently U.S. and Canada only. **Reference promo code: 2018MCTI**

Free Ovation M from choice of the following:

**Mechanical Volume Adjustable
single channel pipette**

Select free	Product #	Volume Range	Color
	1070-0002	0.2 - 2 µL	Red
	1070-0010	1 - 10 µL	Orange
	1070-0020	2 - 20 µL	Yellow
	1070-0100	10 - 100 µL	Green
	1070-0200	20 - 200 µL	Blue
	1070-1000	100 - 1000 µL	Purple



*Alternate option (special for trade-in promo):

I want free Product You can choose to get up to \$250* worth of FREE VistaLab product instead of a free Ovation M pipette.

*Free product value at VistaLab MSRP.
See www.vistalab.com for products and pricing.

Shipping information for FREE pipette

Name

Company/Institution

Address


Phone

Email



Patented ergonomic design for more comfort with superior accuracy & precision

VistaLab Product #		Pipette Color	Volume Range
8 Channel	12 Channel		
1160-0020	1060-0020	Yellow	0.5 – 20µL
1160-0125	1060-0125	Green	2 – 125µL
1160-0250	1060-0250	Blue	5 – 250µL
n/a	1060-0850	Grey	25 – 850µL
1160-1250	n/a	Purple	25 – 1250µL



- **Easy No Pounding Tip Acquisition**

Secure tip acquisition with audible "click-on" confirmation so no pounding required and get equal loading across all channels




- **Easy Low-Force Tip Ejection**

No-thumb-strain, spring-loaded, "flick-off" tip ejection

- **Stands on bench even with tip attached to maintain sterility.**

- **Color coded by volume for easy identification on benchtop**



- **Programmable Functions**

Multi-dispense, mixing, serial dilution, reverse pipetting, and custom volume pre-sets (up to 3)