







VistaLab Product No.	Color	Volume Range	Dispensing Increments	Accuracy	Precision
1070-0002		0.2 – 2 $\mu$ L	0.002 $\mu$ L	$\pm$ 1.5% at 2 $\mu$ L $\pm$ 2.7% at 1 $\mu$ L $\pm$ 9% at 0.2 $\mu$ L	0.9% at 2 $\mu$ L 1.8% at 1 $\mu$ L 7.6% at 0.2 $\mu$ L
1070-0010		1 – 10 $\mu$ L	0.01 $\mu$ L	$\pm$ 1% at 10 $\mu$ L $\pm$ 1.5% at 5 $\mu$ L $\pm$ 2.5% at 1 $\mu$ L	0.4% at 10 $\mu$ L 0.6% at 5 $\mu$ L 1.2% at 1 $\mu$ L
1070-0020		2 – 20 $\mu$ L	0.02 $\mu$ L	$\pm$ 1% at 20 $\mu$ L $\pm$ 1.5% at 10 $\mu$ L $\pm$ 7% at 2 $\mu$ L	0.3% at 20 $\mu$ L 0.5% at 10 $\mu$ L 2% at 2 $\mu$ L
1070-0100		10 – 100 $\mu$ L	0.1 $\mu$ L	$\pm$ 0.8% at 100 $\mu$ L $\pm$ 0.8% at 50 $\mu$ L $\pm$ 3.5% at 10 $\mu$ L	0.15% at 100 $\mu$ L 0.24% at 50 $\mu$ L 1% at 10 $\mu$ L
1070-0200		20 – 200 $\mu$ L	0.2 $\mu$ L	$\pm$ 0.8% at 200 $\mu$ L $\pm$ 0.8% at 100 $\mu$ L $\pm$ 2.5% at 20 $\mu$ L	0.15% at 200 $\mu$ L 0.25% at 100 $\mu$ L 1% at 20 $\mu$ L
1070-1000		100 – 1000 $\mu$ L	1.0 $\mu$ L	$\pm$ 0.8% at 1000 $\mu$ L $\pm$ 0.8% at 500 $\mu$ L $\pm$ 3% at 100 $\mu$ L	0.15% at 1000 $\mu$ L 0.2% at 500 $\mu$ L 0.6% at 100 $\mu$ L
1070-0998	Ovation Mechanical 3-piece Set (Includes 1-10 $\mu$ L, 10-100 $\mu$ L, 100-1000 $\mu$ L)				
1070-0999	Ovation Mechanical 4-piece Set (Includes 0.2-2 $\mu$ L, 2-10 $\mu$ L, 20-200 $\mu$ L, one 100-1000 $\mu$ L) Sets include 3 shelf hangers and 3 racks of tips.)				



Ovation M