	VistaLab Product No.	Color	Volume Range	Dispensing Increments	Accuracy	Precision	
Quick-Set Adjustable Volume Models	1057-0002		0.2 – 2μL	0.002μL	±1.5% at 2μL ±2.7% at 1μL ±9% at 0.2μL	0.9% at 2μL 1.8% at 1μL 7.6% at 0.2μL	Ovation QS
	1057-0010		1 – 10μL	0.01μL	±1% at 10μL ±1.5% at 5μL ±2.5% at 1μL	0.4% at 10μL 0.6% at 5μL 1.2% at 1μL	
	1057-0020		2 – 20μL	0.02μL	±1% at 20μL ±1.5% at 10μL ±7% at 2μL	0.3% at 20μL 0.5% at 10μL 2% at 2μL	
	1057-0100	•	10 – 100μL	0.1μL	±0.8% at 100μL ±0.8% at 50μL ±3.5% at 10μL	0.15% at 100μL 0.24% at 50μL 1% at 10μL	
	1057-0200		20 – 200μL	0.2μL	±0.8% at 200μL ±0.8% at 100μL ±2.5% at 20μL	0.15% at 200μL 0.25% at 100μL 1% at 20μL	
	1057-1000		100 – 1000μL	1.0μL	±0.8% at 1000μL ±0.8% at 500μL ±3% at 100μL	0.15% at 1000μL 0.2% at 500μL 0.6% at 100μL	
	1057-0999 Ovation Quick-Set 3-piece Set (Includes 2-20μL, 20-200μL, 100-1000μL, 3 shelf hangers and 3 racks of tips.)						