

# Pipetting is *FASTER* with **ali-Q<sup>TM</sup> 2**

*Easy aliquoting. More accurate and reproducible than by eye.*

- Faster and easier multidispensing with NO NEED to watch the lines.
- Highly reproducible accuracy with the press of one button!
- Repeat pipetting using any serological pipet.
- Plus conventional pipetting too. ali-Q is a 2-in-1 pipet controller!



ali-Q 2 VS video



*Ideal for cell culture, dispensing media, filling plates, tubes, cryovials, etc.*



**VistaLab<sup>TM</sup>**  
TECHNOLOGIES  
Think Differently. Feel Better.

*"Great for dispensing the same amount repeatedly...You can concentrate on the tip of the pipette and the culture plates because you don't have to watch the graduations on the pipette." - For more reviews see [www.selectscience.net](http://www.selectscience.net)*



# Faster and Easier Aliquoting. No Need to Watch the Lines. More Accurate and Reproducible than by Eye.

## 2-in-1 Repeating plus Conventional Pipetting with a single 2-in-1 Controller

- Push one button to aliquot / multidispense with unmatched repeatability.
- Compatible with all brands and sizes of serological pipets; optimal performance when using Wobble-not™ pipets.
- Innovative system with computer controlled sensors and valves for maximum accuracy and precision.\*  
\*get 2% accuracy and 2% CV at 0.5mL when using Wobble-not™ serological pipets



## ali-Q 2 Series Additional Features

- Variable Speed (VS) model available to accommodate your multidispensing needs.
- On-board Calibration Software
- Removable kickstand
- Lighter weight

Description	Product #	Volume Range
ali-Q 2	2100-0005	0.5 - 5.0 mL
ali-Q 2 LS	2100-0003	0.3 - 3.0 mL
ali-Q 2 VS	2100-1005	0.5 - 5.0 mL

