


Repeat Pipetting is *EASIER* with

 **ali-Q 2**TM
aliquoting pipet controller



• aliquot volume dial


New!
Variable aliquot
Speed (VS) model
now available



• aliquot speed dial



• aliquot button



The ali-Q 2 Series:

Variable aliquoting speed model

- ali-Q 2 VS

Fixed aliquoting speed models

- ali-Q 2 (high speed)
- ali-Q 2 LS (Low Speed)

VistaLab Catalog #	max. aliquot volume
2100-1005	5 mL
2100-0005	5 mL
2100-0003	3 mL

- **Repeat pipetting** / multidispensing with less fatigue
- Highly accurate and reproducible aliquots with one touch
- Use Serological Pipets but no need to watch the lines
- In-lab calibration with on-boarded software