




- ✓ RNase/DNase free
- ✓ Non-pyrogenic
- ✓ No animal by-products
- ✓ Low 'dead' space



## Reagent Reservoirs & Reservoir Coolers: Ideal pipetting companions

We never cut corners, but when re-designing our reagent reservoirs, it was the right thing to do – along with a host of other improvements:

- Trimmed base corners and distinct pour-off spouts make it easy to empty reservoir without spills
- Unique deep well “trough within a trough” design allows even small amounts of liquid to be accessible to all pipette tips, with minimal dead volume
- Sterile, non-pyrogenic, certified RNase/DNase free, and animal free – in bulk or individually wrapped packaging
- Well-marked graduated increments on the inside wall minimize reagent waste
- 12 Channel Reservoir with lid available. Autoclavable.
- Optional lid available for 100 mL units so reagents are protected from evaporation and dust
- Wide base and sturdy construction provides stability and avoids inadvertent spills
- Molded corner structures allow bulk packaged units to be easily separated for use
- Clearly marked for recycling



Here's an innovative way to keep your work from heating up – our reusable Reservoir Cooler:

- Maintains reagents at <math><5^{\circ}\text{C}</math> for up to 2 hours
- Preserve integrity of reagents and extend working time with full reservoir cooling
- Custom fit for the 100 mL reagent reservoir
- Non-toxic blue gel is encased in heavy-duty PVC, chills or freezes quickly, and can be used over and over


# Ordering Information

VistaLab Reagent Reservoirs are available in a variety of volume sizes and packaging configurations.

All VistaLab Reagent Reservoirs are made from Animal-Free raw materials



RESERVOIR SIZE:		MINI	SMALL	SMALL	MEDIUM	LARGE
DESCRIPTIONS	Bulk, 100/case Sterile, non-pyrogenic, RNase/DNase free	10mL <b>3054-1012</b>	25mL <b>3054-1002</b>	25mL Divided <b>3054-1004</b>	50mL <b>3054-1006</b>	100mL <b>3054-1000</b>
	Bulk, 200/case Sterile, non-pyrogenic, RNase/DNase free	10mL <b>3054-2012</b>	25mL <b>3054-2002</b>		50mL <b>3054-2006</b>	100mL <b>3054-2000</b>
	Bulk, 100/case non-sterile	10mL <b>4054-1012</b>	25mL <b>4054-1002</b>	25mL Divided <b>4054-1004</b>	50mL <b>4054-1006</b>	100mL <b>4054-1000</b>
	Individually wrapped, 50/case Sterile, non-pyrogenic, RNase/DNase free	10mL <b>3054-1013</b>	25mL <b>3054-1003</b>	25mL Divided <b>3054-1005</b>	50mL <b>3054-1007</b>	100mL <b>3054-1001</b>
	Individually wrapped, 100/case Sterile, non-pyrogenic, RNase/DNase free	10mL <b>3054-2013</b>	25mL <b>3054-2003</b>		50mL <b>3054-2007</b>	100mL <b>3054-2001</b>
	Individually wrapped with lid, 50/case Sterile, non-pyrogenic, RNase/DNase free				100mL <b>3054-3001</b>	
	Bulk with lid, 100/case Sterile, non-pyrogenic, RNase/DNase free				100mL <b>3054-3000</b>	
	Reservoir Cooler for 100mL reservoir				<b>3054-5000</b>	



10mL Reagent Reservoir



25mL Reagent Reservoir



25mL Divided Reagent Reservoir



50mL Reagent Reservoir

### 12 Channel Reservoirs (polypropylene, autoclavable)

- 3054-1010 12 Channel Reservoir with lid, sterile, individually wrapped, 10/case
- 3054-1011 12 Channel Reservoir with lid, sterile, individually wrapped, 25/case
- 4054-1010 12 Channel Reservoir with lid, non-sterile, bulk, 10/case
- 4054-1011 12 Channel Reservoir with lid, non-sterile, bulk, 25/case



12 Channel Reservoir with Lid (Autoclavable)



100mL Reagent Reservoir Cooler



100mL Reagent Reservoir



100mL Reagent Reservoir with Lid