



# Micro D-Tipper

VistaLab Catalog #	Color	Volume	Accuracy	Precision
1802	Red	2 $\mu$ L	$\pm$ 2.00%	2.50%
1803	Gray	3 $\mu$ L	$\pm$ 1.50%	2.00%
1804	White	4 $\mu$ L	$\pm$ 1.50%	2.00%
1805	Silver	5 $\mu$ L	$\pm$ 1.00%	1.50%
1806	Green	6 $\mu$ L	$\pm$ 1.00%	1.50%
1807	Purple	7 $\mu$ L	$\pm$ 1.00%	1.50%
1808	Black	8 $\mu$ L	$\pm$ 1.00%	1.50%
1809	Blue	9 $\mu$ L	$\pm$ 1.00%	1.50%
1810	Orange	10 $\mu$ L	$\pm$ 1.00%	1.00%

