



The ali-Q pipet controller is an innovative all-in-one product that allows for both conventional and repeat pipetting using serological pipets. No other controller has ali-Q's repeat pipetting capability that eliminates the need to watch the lines on the serological pipet for faster and easier aliquoting than having to aliquot "by eye". Because of ali-Q's computer controlled technology, you get accurate volume dispenses every time, meaning consistent repeatability!

The only 2-in-1 Pipettor

Conventional Pipetting



NEW! for ali-Q 2



3 models

Fixed High Speed
Fixed Low Speed
Variable Speed



LED & Speed Dial

2 multifunction LEDs
Speed dial (for VS model only)



Onboard Calibration Software

Invert to initiate protocol



Removable Kickstand