

Wobble-not™

Low Insertion Force Serological Pipet

**Unique, patented
No Drip Design!**



Do any of your other serological pipets look like this?

2-tier step design

- = 2 zones of contact in your pipet controller = Less force to insert
- = Less "wobble" = Max stability
- = Better seal = Less dripping = Better accuracy

LOSE THE WOBBLE

**LOSE
THE
DRIP**



Ordering Information

Pipet Volume	VistaLab Catalog #	Quantity per case
--------------	--------------------	-------------------

Standard length

Individually wrapped (paper/plastic), Sterile

5mL	4090-0005	200
10mL	4090-0010	200
25mL	4090-0025	200
50mL	4090-0050	100
100mL	4090-0100	50

Individually wrapped (plastic/plastic), Sterile

5mL	4190-0005	200
10mL	4190-0010	200
25mL	4190-0025	200
50mL	4190-0050	100
100mL	4190-0100	50

Bulk packaged (resealable plastic bag), Sterile

5mL	4290-0005	500
10mL	4290-0010	400
25mL	4290-0025	150
50mL	4290-0050	100
100mL	4290-0100	60

Shorty length

Individually wrapped (paper/plastic), Sterile

10mL	5090-0010	200
25mL	5090-0025	100

Why it works

