

CERTIFICATE ID: 227610

PIPETTE IDENTIFICATION

CUSTOMER ID	EQ000214	CATALOG NO.	1057-0002
CUSTOMER NAME	Matthew Smith	MODEL	Ovation QS 0.2-2uL
ADDRESS 1	Accelerate Diagnostics	SERIAL NO.	127768
ADDRESS 2	3950 S. Country Club Road	REPAIR NO.	WD29381
CITY, STATE, ZIP	Tucson AZ 85714	TECHNICIAN ID	Y.B.

TESTING CONDITIONS

WATER TEMP (°C)	20.88	BALANCE MODEL	Sartorius MPS 6.6S
AIR TEMP (°C)	20.50	BALANCE SERIAL NO.	36602317
RELATIVE HUMIDITY (%)	51.7	TIPS USED	4060-1002
BAROMETRIC PRESSURE (mmHg)	742.70	"Z" FACTOR	1.0030
		TEST METHOD	ISO 8655

TESTING INFORMATION

TESTING DATA

(measured weight - mg)	<u>Volume 1</u>	<u>Volume 2</u>	<u>Volume 3</u>
Weighing 1	0.200	1.003	2.003
Weighing 2	0.204	0.999	2.004
Weighing 3	0.203	1.004	2.001
Weighing 4	0.204	1.002	2.003
Weighing 5	0.202	1.002	2.005
Weighing 6	0.204	1.005	2.000
Weighing 7	0.196	1.002	2.007
Weighing 8	0.201	1.005	2.003
Weighing 9	0.204	1.004	2.010
Weighing 10	0.204	1.006	2.002
Mean (mg)	0.202	1.003	2.004

CALIBRATION VERIFICATION RESULTS

Mean Test Volume (µL)	0.203	1.006	2.010
Accuracy (±%)	1.50	0.60	0.50
Precision (CV%)	1.49	0.20	0.15

PERFORMANCE SPECIFICATIONS

Nominal Test Volume (µL)	0.20	1.00	2.00
Accuracy (±%)	9.00	2.70	1.50
Precision (CV%)	7.60	1.80	0.90

PIPETTE STATUS **PASS** **PASS** **PASS**

Comments:

Tested & Authorized By: Yeha Balal

Title: Technician

Date of Testing & Authorization 12/22/2020