

# Ovation<sup>®</sup> M

The Better Way To Pipette

- Fully Mechanical operation
- Comfortable, ergonomic design
- Fully autoclavable




Performance, ergonomics and value  
in one elegant design.

Actual Size



Simple conventional  
volume setting



Easy tip acquisition  
EVERY TIME




Ergonomically-correct posture  
that ensures pipetting comfort

# Innovation you can feel



## Fits your hand like a glove.

Ovation's contoured shape and adjustable hook means a custom fit for every user in the lab.



## Easy tip acquisition & discard.

An audible "click" indicates tips are securely seated and ready for use. And discarding tips takes just the gentle push of a button.



## Light thumb action.

A very low-force plunger allows the user to manually control the speed of aspirating and dispensing with minimal effort.



## Fully autoclavable.

No disassembly needed.




6 models to cover 0.2–1000µL



## Conveniently open PCR and microtubes.

Ovation has a cap opener built right into its base.



## Simple in-lab calibration.

No need to send the pipette out.

## Order Information

Celltreat Catalog #	Color	Volume Range	Dispensing Increments	Accuracy*	Precision*
1070-0002	Red	0.2 – 2µL	0.002µL	±1.5% at 2µL ±2.7% at 1µL ±12% at 0.2µL	0.9% at 2µL 1.8% at 1µL 7.6% at 0.2µL
1070-0010	Orange	1 – 10µL	0.01µL	±1% at 10µL ±1.5% at 65µL ±2.5% at 2µL	0.4% at 10µL 0.6% at 5µL 1.2% at 1µL
1070-0020	Yellow	2 – 20µL	0.02µL	±1% at 20µL ±1.5% at 10µL ±7% at 2µL	0.3% at 20µL 0.5% at 10µL 2% at 2µL
1070-0100	Green	10 – 100µL	.1µL	±0.8% at 100µL ±0.8% at 50µL ±3.5% at 10µL	0.15% at 100µL 0.24% at 50µL 1% at 10µL
1070-0200	Blue	20 – 200µL	.2µL	±0.8% at 200µL ±0.8% at 100µL ±2.5% at 20µL	0.15% at 200µL 0.25% at 100µL 1% at 20µL
1070-1000	Purple	100 – 1000µL	1.0µL	±0.8% at 1000µL ±0.8% at 500µL ±3% at 100µL	0.15% at 1000µL 0.2% at 500µL 0.6% at 100µL

\* Specifications subject to change



## VistaLab tips – A fitting companion.

Whether you require tips that are sterile, certified RNase/ DNase free, non-pyrogenic, trace metal certified or filtered, there's a genuine VistaLab tip to meet your needs. Visit [www.celltreat.com](http://www.celltreat.com) for complete information on all the tip options for Ovation pipettes.

## Ovation M multi-pipette sets

3 and 4 multi-pipette sets are available at special pricing.

Cat# 1070-0998 - 1-10µL, 10-100µL and 100-1000µL

Cat# 1070-0999 - 0.2-2µL, 2-20µL and 20-200µL and 100-1000µL



## VistaLab Technologies, Inc.

2 Geneva Road (888) 652-6520  
Brewster, NY 10509 USA (914) 244-6226

©2015 VistaLab Technologies, Inc. All rights reserved.

Printed in the U.S.A. Lit. #374062 Rev D