

Pipette Tips

Find your pipette model and size in chart and appropriate tip size at the top of the column

	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL	5 mL Tips	10 mL Tips
Ovation ESC, E8 & E12	20 µL	125 µL 250 µL	850 µL 1250 µL		
Ovation QS & M	2 µL 10 µL 20 µL	100 µL 200 µL	1000 µL		
Ovation F1 & F2	1 to 20 µL	21 to 200 µL	201 to 1000 µL		
Ovation Macro5 & Macro10				5 mL	10 mL
All MLA models		2 to 200 µL	201 to 1000 µL	1 mL to 5 mL	6 mL to 10 mL

Non-Sterile Tips

Description	Celltreat Catalog #		
	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL
VistaRak™ 192 tips/rack, 5 racks	4060-1002	4060-2004	4060-3004
VistaRak, Low Retention 192 tips/rack, 5 racks	4070-1002LR	4070-2004LR	4070-3004LR
VistaStak™ 192 tips/layer 5 layers for Small tips (960 tips) or 3 layers for Micro & Large tips (576 tips)	4060-9024	4060-9025	4060-9026
Stacked Rack 200 tips/layer 5 layers for Small tips 3 layers for Large tips		9025 ¹	9026 ¹
Stacked Rack Trace Metal Certified 200 tips/layer 5 layers for Small tips 3 layers for Large tips		9022 ¹	9023 ¹
VistaBulk™ 1000 tips/bag	4058-1000	4058-2000	4058-3000
Protectainer™ Bulk Pack 1000 Small tips or 750 Large tips		4025	4026
Econo-Pak™ Bulk Pack 1000 tips		4225	4226

Note: 1. Not for use with Ovation Multichannel pipettes.

Sterile

Description	Celltreat Catalog #		
	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL
VistaRak, Sterile Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4060-1032	4060-2132	4060-3132
VistaRak, Sterile, Low Retention Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4070-1032LR	4070-2132LR	4070-3132LR
VistaTip™, Individually wrapped, Sterile 200 tips		2025	2026
VistaTip, Individually wrapped, Sterile Certified RNase/DNase free, non-pyrogenic 200 tips		2027	2028

Sterile Filtered

Description	Celltreat Catalog #		
	Micro Tips 25 µL	Small Tips 250 µL	Large Tips 1250 µL
VistaRak, Sterile, Filtered Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4060-1332 ² 4060-1333 ¹	4060-2332	4060-3332
VistaRak, Sterile, Filtered, Low Retention Certified RNase/DNase free, non-pyrogenic 192 tips/rack, 5 racks	4070-1332LR ² 4070-1333LR ¹	4070-2332LR	4070-3332LR

Notes: 2. Not for use with Ovation QS or Ovation M 2 uL or Ovation F2 1 uL or 2 uL.

¹ 2 uL tips for Ovation 1 uL and 2 uL pipettes

All Macro Volume Tips

Description	Celltreat Catalog #	
	Macro Tips 5 mL*	Macro Tips 10 mL
Box, Non-Sterile 100 tips/box for 5 mL 35 tips/box for 10 mL	9048	9050
VistaClear™ Box, Non-Sterile 35 tips/box		4058-6100
VistaBulk, Non-Sterile 250 tips/bag for 5 mL 100 tips/bag for 10 mL	4058-5000	4058-6000
VistaClear Box, Sterile Certified RNase/DNase free, non-pyrogenic 60 tips/box for 5 mL 35 tips/box for 10 mL	4058-5102	4058-6102
VistaClear Box, Sterile, Filtered Certified RNase/DNase free, non-pyrogenic 60 tips/box for 5 mL 35 tips/box for 10 mL	4058-5133	4058-6133
VistaTip, Individually wrapped, Sterile, Filtered Certified RNase/DNase free, non-pyrogenic 50 tips	4058-5332	4058-6332

*All 5 mL macro tips are graduated at 0.1 mL increments.