

Catalog #	Volume	Accuracy	Precision
2100-0003	3.0 mL	±2.0%	±1.0%
	0.3 mL	±3.0%	±3.0%
2100-0005	5.0 mL	±2.0%	±1.0%
	0.5 mL	±2.0%	±2.0%
2100-1005*	5.0 mL	±2.0%	±1.0%
	0.5 mL	±2.0%	±2.0%
2100-1010*	10.0 mL	±2.0%	±1.0%
	1.0 mL	±2.0%	±2.0%

*Performed on speed 7

**Performance specs based on use with Wobble-not brand serological pipets. Performance may vary slightly when used with other brands of serological pipets.

