

Reagent Reservoirs

- Proprietary surface treatment for lowest residual reagent/highest reagent recovery
- Trough-within-a-trough design to recover the “last few drops”
- Pour-off spouts at every corner
- A wider base for stability
- Free of animal products

Description	Celltreat Catalog #				
	10 mL	25 mL	25 mL Divided	50 mL	100 mL
Bulk, 100/case Non-sterile	4054-1012	4054-1002	4054-1004	4054-1006	4054-1000
Bulk, 100/case Sterile, non-pyrogenic, RNase/DNase free	3054-1012	3054-1002	3054-1004	3054-1006	3054-1000
Bulk, 200/case Sterile, non-pyrogenic, RNase/DNase free	3054-2012	3054-2002		3054-2006	3054-2000
Individually wrapped, 50/case Sterile, non-pyrogenic, RNase/DNase free	3054-1013	3054-1003	3054-1005	3054-1007	3054-1001
Individually wrapped, 100/case Sterile, non-pyrogenic, RNase/DNase free	3054-2013	3054-2003		3054-2007	3054-2001
Individually wrapped with lid, 50/case Sterile, non-pyrogenic, RNase/DNase free					3054-3001
Lid for reservoir, 25/bag, 100/case Sterile, non-pyrogenic, RNase/DNase free					3054-3000
Reservoir Cooler for 100 mL reservoir					3054-5000

12-channel Reservoirs

- with Lid
- Polypropylene
- Autoclavable

Bulk, 10/case, Non-sterile
 Bulk, 25/case, Non-sterile
 Individually wrapped, 10/case, Sterile
 Individually wrapped, 25/case, Sterile

4054-1010
 4054-1011
 3054-1010
 3054-1011

